

TurboReader™



Veterinary
Point-of-Care system



Diagnostic point-of-care reader which together with single-use reagent cuvettes performs one analysis within 3 minutes.

www.euris.org


Turbo Reader™

Veterinary
Point-of-Care system




The TurboReader™ is an easy-to-use point-of-care system for rapid veterinary diagnostics. It provides the accuracy of the large-scale laboratories at patient care level in veterinary hospitals and during field use.

An expanding range of tests (for both serum/plasma and whole blood) are in development. The Canine C-Reactive Protein (CRP) TurboReader™ Assay is intended for *in vitro* determination of CRP in dogs. Canine CRP is a real-time marker for systemic inflammation, as well as, response to ongoing clinical treatment. Its normal plasma concentration in healthy dogs is <35 mg/L.



The rapid measuring procedure involves the pipetting of 20 µl of sample into a cuvette containing reaction buffer R1. After 1 minute, 3 drops of reaction buffer R2 are added. The quantitative result is presented within 3 minutes on the touch screen color display. High accuracy and sensitivity are obtained.

Excellent correlation with validated method: The assay performance has been compared with a validated reference method on 19 dog serum samples.



Connecting of an optional barcode reader enables patient traceability. Results can be printed with a thermal printer. Connection to PC or LIM system will be available and requires a software upgrade.

Order information

Product	Article No.
TurboReader™ Instrument	2200-04
Canine CRP TurboReader™ Assay	2530-01

European Institute of Science AB

Scheelevägen 27, 4th Floor SE-223 70 Lund, Sweden

Phone: +46 (46) 286 2230

info@euris.org | www.euris.org

