

# Equine Serum Amyloid A (SAA) TurboReader™ Assay



GET THE CLEAR  
CLINICAL PICTURE  
WITH A SIMPLE TEST:

---

- ✓ **DETECT PRESENCE &  
QUANTIFY SEVERITY  
OF INFECTION**
  - ✓ **MONITOR RESPONSE  
TO ANTIBIOTICS  
& SURGERY**
  - ✓ **SCREENING TOOL  
FOR INFECTIONS**
- 





# DETECTION AND REAL-TIME MONITORING OF SYSTEMIC INFLAMMATION IN HORSES



## Equine SAA TurboReader™ Assay

Article No. 2534-01

Sample type: ..... Serum or Li-heparin Plasma  
Sample volume: ..... 5 µl  
Measuring range: ..... 10 – 600 mg/l  
Kit size: ..... 20 tests  
Analysis time: ..... 9.5 minutes  
Clinical decision limit: ..... 30 mg/l

**\*Horse specific calibrated assay.**